

Oligonucleotide research has exploded in recent years, requiring new and sophisticated techniques not only in their application, but also analysis. Phosphorothioate modification of oligonucleotides (S-oligos) is a technique used to improve stability and slow degradation by nucleases. In the synthesis of S-oligos it is imperative that modifications are as directed, whether it be a partially or fully modified S-oligo. [Cosmocore C18](#) columns are shown here in a method to elucidate minute differences in chemical modification with excellent resolution.

Acquisition Method

Column: [Cosmocore C18](#)
Column Size: 2.1 x 100 mm
Mobile Phase: A: 100mM TEAA water; B: Acetonitrile
Gradient: 10 → 25 %B (0 → 10 min)
Flow Rate: 0.4 mL/min
Temperature: 50 Celsius
Detection: UV 260
Sample: Phosphorothioate modified DNA oligos
Inj. Vol.: 1 µL

Sample Sequences

A: 20-mer All PS

T*T

B: 20-mer 1PO 5'

TT*T*T*T*T*T*T*T*T*T*T*T*T*T*T*T*T*T

C: 20-mer 2PO 5'+ middle

TT*T*T*T*T*T*T*T* T*TT*T*T*T*T*T*T*T

T* = phosphorothioate linkage




Figure 1: Overlaid chromatograms of fully and partially phosphorothioated DNA oligos. Excellent resolution is observed between oligos of the same length and differing number of modifications.