

HPLC column

ChromaNik
ChromaNik Technologies Inc.

SunBridge




ChromaNik Technologies

Nacalai USA
Innovations for Life Sciences.



SunBridge C18

Ethylene cross-linked silica gel




Synthesis of hybrid silica gel


Particle size: 5 μm , Pore volume: 0.8 mL/g,
Specific surface area: 190 m^2/g , Pore diameter: 15 nm

Bonding with Trifunctional C18

Double End-capping




Elevated temperature end-capping



Hybrid silica gel + special end-capping

Evaluation of hydrogen bonding, hydrophobicity and steric selectivity



As a comparison, Hybrid C18 manufactured by W Company, which is ethylene crosslinked silica gel C18 in which an ethylene chain is incorporated into the silica skeleton, was used.

Column: Company W Hybrid C18 3.5 μm , 5 μm ,

SunBridge C18 5 μm ,



Column dimensions: 150 x 4.6 mm

Mobile phase: $\text{CH}_3\text{OH}/\text{H}_2\text{O}=75/25$

Flow rate: 1.0 mL/min

Temperature: 40 $^\circ\text{C}$


Sample: 1 = Uracil, 2 = Caffeine, 3 = Phenol, 4 = Butylbenzene,



5 = o-Terphenyl,



6 = Amylbenzene,




7 = Triphenylene



	Hydrogen bonding $\alpha(\text{Caffeine/Phenol})$	Hydrophobicity $\alpha(\text{Amylbenzene/Butylbenzene})$	Steric selectivity $\alpha(\text{Triphenylene/o-Terphenyl})$	Specific surface area (m^2/g)	Carbon loading (%)
Company W Hybrid C18 3.5 μm	0.38	1.54	1.33	185	17.4%
Company W Hybrid C18 5 μm	0.39	1.52	1.33	188	18.0%
SunBridge C18 5 μm	0.42	1.54	1.32	190	16.1%

Comparison of peaks for pyridine, oxine, and nortriptyline



Column dimension:

150 x 4.6 mm for 3.5 and 5 μm
100 x 2.1 mm for 1.7 μm

Mobile phase: CH₃OH/H₂O=30/70

Flow rate: 1.0 mL/min for 3.5 and 5 μm
0.2 mL/min for 1.7 μm

Temperature: 40 °C

Detection: UV@250 nm

Sample: 1 = Uracil

2 = Pyridine

3 = Phenol



Column dimension:

150 x 4.6 mm for 3.5 and 5 μm
100 x 2.1 mm for 1.7 μm

Mobile phase: CH₃CN/20mM H₃PO₄=10/90

Flow rate: 1.0 mL/min for 3.5 and 5 μm
0.2 mL/min for 1.7 μm

Temperature: 40 °C

Detection: UV@250 nm

Sample: 1 = 8-Quinololin (Oxine)

2 = Caffeine



Column dimension:

150 x 4.6 mm for 3.5 and 5 μm
100 x 2.1 mm for 1.7 μm

Mobile phase: Acetonitrile/10mM ammonium acetate pH6.8=40/60



Flow rate: 1.0 mL/min for 3.5 and 5 μm
0.3 mL/min for 1.7 μm

Temperature: 40 °C

Detection: UV@250 nm

Sample: 1=Uracil, 2=Propranolol,


3= Nortriptyline, 4=Amitriptyline



Company W Hybrid C18 showed severe tailing of metal-chelating, oxine.

Compared to Company W Hybrid C18, SunBridge C18 showed a higher number of plates and less tailing of peaks for metal-chelating and basic compounds.

Comparison of Stability under pH 11.5 and pH 1.0 conditions




Durable test condition (pH 11.5)

Column dimension: 150 x 4.6 mm
Mobile phase: Methanol/50mM potassium phosphate pH 11.5=10/90
Flow rate: 1.0 mL/min
Temperature: 40 °C

Measurement condition

Mobile phase: Acetonitrile/water=70/30
Flow rate: 1.0 mL/min
Temperature: 40 °C
Sample: 1 = Uracil, 2 = Butylbenzene



Durable test condition (pH 1.0)

Column dimension: 50 x 2.1 mm
Mobile phase: Acetonitrile/1.0% TFA pH1=10/90
Flow rate: 0.4 mL/min
Temperature: 80 °C

Measurement condition

Same as above

SunBridge C18 showed almost the same or higher stability under pH 11.5 condition as Company W Hybrid C18, which shows an order of magnitude higher stability than other companies' hybrid C18. In addition, the decrease in retention due to alkali was suppressed by end-capping with a bidentate silylating reagent.

Ordering information

Packings	Inner diameter (mm)	2.1	3.0	4.6	10	20	US P category
	Length (mm)	Catalog number	Catalog number	Catalog number	Catalog number	Catalog number	
SunBridge C18, 5 μm	50	JB3941	JB3341	JB3441	-----	-----	L1
	100	JB3961	JB3361	JB3461	-----	-----	
	150	JB3971	JB3371	JB3471	-----	-----	
	250	JB3981	JB3381	JB3481	JB3781	JB3881	

* SunBridge columns packed with 3 μm and 1.8 μm particles will be released sequentially.

Nacalai USA, Inc.

6625 Top Gun Street, Suite 107, San Diego, CA 92121
 TEL: 858 404 0403
 E-mail: info@nacalaiusa.com
 URL: www.nacalaiusa.com

ChromaNik Technologies Inc.

6-3-1 Namiyoke, Minato-ku, Osaka, 552-0001 Japan
 TEL: +81-6-6581-0885 FAX: +81-6-6581-0890
 E-mail: info@chromanik.co.jp
 URL: www.chromanik.co.jp